

Product Datasheet

Normalized Human Serum H1091

Product Name	Normalized Human Serum
Catalog Number	H1091
Description	Normalized Human Serum
Tested Applications	QC Materials, Base Matrix
Alternative Names	Plasma Derived Serum, Defibrinated Plasma, Filtered Plasma
Brand	Bioresource Technology
Form/Appearance	Liquid
Storage	-20°C
Note	This product is final formulated to the following parameters: Total Protein: Reported in g/dL Sodium: Reported in mmol/L Potassium: Reported in mmol/dL Chloride: Reported in mmol/dL Cholesterol: Reported in mg/dL Triglycerides: Reported in mg/dL Phosphorus: Reported in mg/dL Creatinine: Reported in mg/dL Calcium: Reported in mg/dL Albumin: Reported in g/dL pH: Reported Microbial Count: <100 CFU/mL
Purity	0.2 micron filtered plasma
Source	Pooled Human Donors