

Product Datasheet

BioCon II - Defibrinated Plasma H2000

Product Name	BioCon II - Defibrinated Plasma
Catalog Number	H2000
Description	BioCon II - Defibrinated Plasma
Tested Applications	QC Materials, Base Matrix
Alternative Names	Plasma Derived Serum, Filtered Plasma, Processed Plasma
Brand	Bioresource Technology
Form/Appearance	Liquid
Storage	-20°C
Note	This product is final formulated to the following parameters: Total Protein: >5.0 g/dL Sodium: 80-150 mmol/L Potassium: 2.0-6.0 mmol/L Chloride: 65-120 mmol/L Cholesterol: <200 mg/dL Triglycerides: <200 mg/dL pH: 7.0-7.6 Bioburden: <100 cfu/mL
Purity	0.2 micron filtered plasma
Source	Pooled Human Donors