

Product Datasheet

Human Control Base Matrix H2001

Product Name	Human Control Base Matrix
Catalog Number	H2001
Description	Human Control Base Matrix
Tested Applications	QC Materials, Base Matrix
Alternative Names	Processed Plasma, Defibrinated Plasma, Filtered Plasma
Brand	Bioresource Technology
Form/Appearance	Liquid
Storage	-20°C
Note	Total Protein: Reported in g/dL Sodium: Reported in mmol/L Potassium: Reported in mmol/L Chloride: Reported in mmol/L Cholesterol: Reported in mg/dL Triglycerides: Reported in mg/dL Phosphorus: Reported in mg/dL Creatinine: Reported in mg/dL Glucose: Reported in mg/dL Calcium: Reported in mg/dL Albumin: Reported in g/dL pH: Reported Microbial Count: Reported in cfu/mL
Purity	0.2 micron filtered plasma
Source	Pooled Human Donors