

## Product Datasheet

### Human Control Base Matrix H2001

<b>Product Name</b>	Human Control Base Matrix
<b>Catalog Number</b>	H2001
<b>Description</b>	Human Control Base Matrix
<b>Tested Applications</b>	QC Materials, Base Matrix
<b>Alternative Names</b>	Processed Plasma, Defibrinated Plasma, Filtered Plasma
<b>Brand</b>	Bioresource Technology
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20°C
<b>Note</b>	Total Protein: Reported in g/dL Sodium: Reported in mmol/L Potassium: Reported in mmol/L Chloride: Reported in mmol/L Cholesterol: Reported in mg/dL Triglycerides: Reported in mg/dL Phosphorus: Reported in mg/dL Creatinine: Reported in mg/dL Glucose: Reported in mg/dL Calcium: Reported in mg/dL Albumin: Reported in g/dL pH: Reported Microbial Count: Reported in cfu/mL
<b>Purity</b>	0.2 micron filtered plasma
<b>Source</b>	Pooled Human Donors