

Product Datasheet

Universal Negative Base Matrix H2420

| | |
|----------------------------|---|
| Product Name | Universal Negative Base Matrix |
| Catalog Number | H2420 |
| Description | Universal Negative Base Matrix |
| Tested Applications | QC Materials, Base Matrix |
| Alternative Names | Processed Plasma, Defibrinated Plasma, Filtered Plasma |
| Brand | Bioresource Technology |
| Form/Appearance | Liquid |
| Storage | -20°C |
| Note | This product is final formulated to the following parameters: Total Protein: 6.0-8.3 g/dL Sodium: <150 mmol/L aHBs: Negative pH: 7.0-7.6 Bioburden: <100 cfu/mL |
| Purity | 0.2 micron filtered plasma |
| Source | Pooled Human Donors |