

## Product Datasheet

### BioCon II CD - Charcoal Stripped Plasma H2550

<b>Product Name</b>	BioCon II CD - Charcoal Stripped Plasma
<b>Catalog Number</b>	H2550
<b>Description</b>	BioCon II CD - Charcoal Stripped Plasma
<b>Tested Applications</b>	QC Materials, Base Matrix
<b>Alternative Names</b>	Processed Plasma, Defibrinated Plasma, Filtered Plasma
<b>Brand</b>	Bioresource Technology
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20°C
<b>Purity</b>	0.2 micron filtered plasma
<b>Source</b>	Pooled Human Donors