

Product Datasheet

Cortisol-3-HRP Antigen 170053

Product Name	Cortisol-3-HRP Antigen
Catalog Number	170053
Description	Synthetic Cortisol
Brand	Medix Biochemica
Form/Appearance	Liquid; Antigen-HRP conjugate
Storage	Storage at 2-8 °C. Do not freeze.
Note	This product was previously available with catalog code P91-92-91H. It has been produced using the same conjugation protocol.
Buffer	Tris base/FBS (1:1), pH 7.0, 0.0016% gentamicin sulfate, 0.005% thimerosal
Matched Pairs	Anti-Cortisol 12901 SPTN-5 100904,
Shelf Life	Unspecified