

Product Datasheet

Anti-h NfL 12601 SPTN-5
100984

Product Name	Anti-h NfL 12601 SPTN-5
Catalog Number	100984
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA, CLIA
Alternative Names	Alzheimer's disease (AD), Traumatic brain injury (TBI), Parkinson's disease (PD), Huntington's disease, Multiple sclerosis (MS). Amyotrophic lateral sclerosis (ALS)
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 95 %
Associated Products	NfL antigen LA666
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
IEF Profile	6.6-7.3
Cross Reactivity	Does not recognize NfM, NfH, NSE or S100B. Recognizes bovine NfL.
Specificity	Antibody recognizes human neurofilament light
Matched Pairs	NfL 12603 and 12604
Shelf Life	Unspecified