

Product Datasheet

Immunoglobulin M (IgM) 340-39

Product Name	Immunoglobulin M (IgM)
Catalog Number	340-39
Alternative Names	IgM Concentrate, Partially Purified IgM
Brand	Lee Biosolutions
Concentration	≥ 15 mg/mL (Dimension® Clinical Chemistry System)
Storage	-20 °C
Purity	Semi-pure
Source	Human Plasma (defibrinated) from Multiple Donors
Buffer	20 mM sodium phosphate, 0.5 M NaCl with 0.05% Proclin 300, pH 7.2,
Shelf Life	1 year