

Product Datasheet

Whole Blood, Pooled Human Donors 991-50-P

Product Name	Whole Blood, Pooled Human Donors
Catalog Number	991-50-P
Brand	Lee Biosolutions
Form/Appearance	Liquid (Red)
Storage	2-8°C
Note	Negative/non-reactive for HIV 1/2, HBsAg, RPR and HCV by FDA approved methods. Inquire for custom anticoagulants and donor pooling.
Source	Pooled Human Donors (≥ 3 donors per pool)
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	K3EDTA as an anticoagulant