

Product Datasheet

Whole Blood, Pooled Donors 991-50-P

Product Name	Whole Blood, Pooled Donors
Catalog Number	991-50-P
Description	Inquire for custom anticoagulants and donor pooling.
Tested Applications	QC Materials
Alternative Names	Blood, whole, total blood, unseparated, complete, fresh, venous
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	2-8°C
Note	Negative/non-reactive for HIV 1/2, HBsAg, RPR and HCV by FDA approved methods. Inquire for custom anticoagulants and donor pooling.
Source	Pooled Human Donors
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	K3EDTA (anticoagulant), inquire for other options