

Product Datasheet

Whole Blood, Single Donor 991-50-S

Product Name	Whole Blood, Single Donor
Catalog Number	991-50-S
Description	Each unit is typically 500mL per donor; however, this varies by unit. If unit volumes are critical, please communicate when ordering. Inquire for custom volumes per donor or custom criteria.
Tested Applications	QC Materials
Alternative Names	Blood, whole, total blood, unseparated, complete, fresh, venous
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	2-8°C
Note	One unit is ~ 500 mL. Smaller volumes are available upon request.
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	K3EDTA (anticoagulant), inquire for other options