

# Product Datasheet

## Whole Blood, Single Human Donor 991-50-S

Product Name	Whole Blood, Single Human Donor
Catalog Number	991-50-S
Description	Donor Age: Reported Donor Race: Reported Donor Gender: Reported
Brand	Lee Biosolutions
Form/Appearance	Liquid (Red)
Storage	2-8°C
Note	One unit is ~ 500 mL. Smaller volumes are available upon request.
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	K3EDTA as an anticoagulant. Inquire for custom anticoagulants.