

Product Datasheet

Pleural Fluid, Single Donor 991-38-S

Product Name	Pleural Fluid, Single Donor
Catalog Number	991-38-S
Description	Donor Age: Reported Donor Race: Reported Donor Gender: Reported
Tested Applications	QC Materials
Alternative Names	Lung, lung cavity, lung cavity fluid, chest, chest fluid, thoracic, pleural effusion, pleural
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	-20°C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None added, inquire if needed