

Product Datasheet

Pleural Fluid, Pooled Human Donors 991-38-P

Product Name	Pleural Fluid, Pooled Human Donors
Catalog Number	991-38-P
Brand	Lee Biosolutions
Form/Appearance	Raw, unprocessed liquid
Storage	-20°C
Source	Pooled Human Donors (≥ 3 donors per pool)
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None