

Product Datasheet

Anti-Human CD117 Monoclonal Antibody, Azide Free Clone B-N20 853.730.005

Product Name Anti-Human CD117 Monoclonal Antibody, Azide Free Clone B-N20

Catalog Number 853.730.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide Free

Tested Applications ELISA

Immunogen Recombinant c-kit

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm.
Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the Stem Cell Factor Receptor, SCF R, c-kit, a 145 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells