

---

## Product Datasheet

### Human Chorionic Gonadotropin 189-12

<b>Product Name</b>	Human Chorionic Gonadotropin
<b>Catalog Number</b>	189-12
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	2-8 °C (short-term) -20°C (long-term)
<b>Purity</b>	≥ 95% (SDS-PAGE)
<b>Source</b>	Human Pregnancy Urine
<b>Molecular Weight</b>	36,700
<b>Biological Activity</b>	≥ 10,000 IU/mg (RID)