

---

## Product Datasheet

### Low Density Lipoprotein (LDL) - MDA Modified 360-27

<b>Product Name</b>	Low Density Lipoprotein (LDL) - MDA Modified
<b>Catalog Number</b>	360-27
<b>Alternative Names</b>	MDA modified LDL
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	≥ 1 mg/mL (Coomassie)
<b>Storage</b>	2-8 °C
<b>Purity</b>	≥ 95% (SDS-PAGE)
<b>Source</b>	Human Plasma
<b>Molecular Weight</b>	2.75 million (ref. Tetz second Edition)