

Product Datasheet

Low Density Lipoprotein (OxLDL) - Copper Oxidized 360-31

Product Name	Low Density Lipoprotein (OxLDL) - Copper Oxidized
Catalog Number	360-31
Alternative Names	OxLDL, Ox-LDL
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	5 - 8 mg/mL (Modified Lowry)
Storage	2-8 °C
Purity	≥ 95% (SDS-PAGE before oxidation)
Source	Human Plasma
Triglyceride Ratio	≤ 0.50 (Triglycerides/Total Cholesterol)