

## Product Datasheet

### High Density Lipoprotein (HDL) 361-25

<b>Product Name</b>	High Density Lipoprotein (HDL)
<b>Catalog Number</b>	361-25
<b>Alternative Names</b>	HDL, HDL-C, HDL Concentrate, HDL Fraction, alpha-Lipoprotein
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1.0 - 5.0 mg/mL (Lowry)
<b>Storage</b>	2-8 °C
<b>Purity</b>	≥ 95% (SDS-PAGE)
<b>Source</b>	Human Plasma
<b>Molecular Weight</b>	175,000-360,000 (Tetz second Edition)