

# Product Datasheet

## Plasma, Single Donor, Normal Recovered K3EDTA 991-58-S-K3ED

Product Name	Plasma, Single Donor, Normal Recovered K3EDTA
Catalog Number	991-58-S-K3ED
Description	Each unit is typically 250-300mL per donor; however, this varies by unit. If unit volumes are critical, please communicate when ordering. Inquire for custom volumes per donor or custom criteria.
Tested Applications	QC Materials
Alternative Names	Blood plasma, plasma, whole plasma, fraction, plasma fraction, anticoagulated
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	-20°C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	Typically EDTA (anticoagulant), inquire for other options