

Product Datasheet

Carbapenemase (OXA-48) Antibody HM1433

Product Name	Carbapenemase (OXA-48) Antibody
Catalog Number	HM1433
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Alternative Names	OXA-48
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM1434
Preservatives	0.09% Sodium Azide