

Product Datasheet

Carbapenemase (KPC) Antibody HM1441

Product Name	Carbapenemase (KPC) Antibody
Catalog Number	HM1441
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Alternative Names	KPC
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM1440 or HM1442
Preservatives	0.09% Sodium Azide