

Product Datasheet

Progesterone-3-BSA Antigen LA814

Product Name	Progesterone-3-BSA Antigen
Catalog Number	LA814
Conjugation	low hapten:BSA ratio (5:1)
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM524, HM525 or HM1335 (detection)
Preservatives	0.01% Sodium Azide