

Product Datasheet

Anti-h Tau R13303 SPTN-5 140040

Product Name	Anti-h Tau R13303 SPTN-5
Catalog Number	140040
Description	Recombinant mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.
Tested Applications	FIA
Alternative Names	Neurology, Mild cognitive impairment (MCI), Alzheimer's disease (AD)
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Epitope	A non-phosphorylated epitope within the overlapping area of the six different CNS Tau variants (352-441 amino acids long)
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.0-6.8
Specificity	Antibody recognizes human total Tau
Shelf Life	Unspecified, storage at 2-8 °C