
Product Datasheet

Recombinant TSH antigen 610175

Product Name	Recombinant TSH antigen
Catalog Number	610175
Species/Host	CHO Cells
Alternative Names	Thyroid Stimulating Hormone (TSH) recombinant human antigen
Brand	Medix Biochemica
Form/Appearance	Liquid
Storage	Store at -20 °C
Source	Recombinant
Buffer	50 mM Tris-HCl, 150 mM NaCl, pH 7.5
Molecular Weight	25 kDa
Expression System	CHO