

## Product Datasheet

### Recombinant TSH antigen 610175

<b>Product Name</b>	Recombinant TSH antigen
<b>Catalog Number</b>	610175
<b>Species/Host</b>	CHO Cells
<b>Alternative Names</b>	Thyroid Stimulating Hormone (TSH) recombinant human antigen
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Lot specific, $\geq 0.5$ mg/ml
<b>Storage</b>	Store at -20 °C
<b>Purity</b>	$\geq 90$ % (HPLC-SEC)
<b>Source</b>	Recombinant
<b>Buffer</b>	50 mM Tris-HCl, 150 mM NaCl, pH 7.5
<b>Molecular Weight</b>	25 kDa
<b>Matched Pairs</b>	Anti-h TSH 5405, 5408, 5409 and Anti-h Alpha subunit 6601
<b>Expression System</b>	Mammalian Cell (CHO)