

Product Datasheet

Anti-Fentanyl 13101 SPTN-5 100998

Product Name	Anti-Fentanyl 13101 SPTN-5
Catalog Number	100998
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/-10 %)
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 95 %
Associated Products	170040
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.6-7.1
Cross Reactivity	Recognizes acetyl fentanyl Does not recognize norfentanyl Risperidone ≤ 1 % 9-hydroxyrisperidone ≤ 1 %
Specificity	Antibody recognizes fentanyl
Matched Pairs	Fentanyl-BSA Antigen 170040 (capture)
Shelf Life	Unspecified