

Product Datasheet

Anti-Human CD44 Monoclonal Antibody, Azide Free Clone B-R8 853.790.005

Product Name Anti-Human CD44 Monoclonal Antibody, Azide Free Clone B-R8

Catalog Number 853.790.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide Free

Tested Applications ELISA

Immunogen U266 cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the Hyaluronate Cell Adhesion Molecule (H-CAM), a 85-250 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells