

Product Datasheet

Serum, Single Donor, Triglyceride 991-24-PS-TRIG

Product Name	Serum, Single Donor, Triglyceride
Catalog Number	991-24-PS-TRIG
Description	Each sample is a minimum of 1mL per donor. Inquire for larger volumes per donor or custom criteria.
Tested Applications	QC Materials
Alternative Names	Serum, fraction, clotted, post-coagulation, coagulated, powder, sera, blood serum, defibrinated, triglyceride
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	991-24-PS-TRIG1 - Low Trig: < 249 mg/dL (Chemiluminescence); 991-24-PS-TRIG2 - Normal Trig: 250-500 mg/dL; 991-24-PS-TRIG3 - High Trig: > 501 mg/dL; Note: Please specify range upon order placement
Storage	-20°C
Note	991-24-PS-TRIG1 - Low Trig: < 249 mg/dL (Chemiluminescence); 991-24-PS-TRIG2 - Normal Trig: 250-500 mg/dL; 991-24-PS-TRIG3 - High Trig: > 501 mg/dL; Note: Please specify range upon order placement
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None added, inquire if needed