

Product Datasheet

Xylazine Antibody HM1435

Product Name	Xylazine Antibody
Catalog Number	HM1435
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Xylazine-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Xylazine 100%, 3-hydroxy xylazine 150%, 4-hydroxy xylazine 180% and 4-hydroxy xylazine-O-glucuronide 200%
Specificity	Xylazine and metabolites
Matched Pairs	LA815
Preservatives	0.05% Sodium Azide