

## Product Datasheet

### Xylazine Antibody HM1437

<b>Product Name</b>	Xylazine Antibody
<b>Catalog Number</b>	HM1437
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF
<b>Immunogen</b>	Xylazine-KLH
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Ion Exchange Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
<b>Cross Reactivity</b>	Xylazine 100%; 3-hydroxy xylazine, 4-hydroxy xylazine and 4-hydroxy xylazine-O-glucuronide 150%
<b>Specificity</b>	Xylazine and metabolites
<b>Matched Pairs</b>	LA815
<b>Preservatives</b>	0.05% Sodium Azide