

Product Datasheet

Xylazine-BSA Antigen LA815

Product Name	Xylazine-BSA Antigen
Catalog Number	LA815
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	≥90% by SDS-PAGE
Buffer	10 mM Phosphate Buffer, pH 7.2
Matched Pairs	HM1435, HM1436, HM1437
Preservatives	0.01% Sodium Azide