

Product Datasheet

Anti-EtG R12403 SPTNZ-1 140052

| | |
|----------------------------|---|
| Product Name | Anti-EtG R12403 SPTNZ-1 |
| Catalog Number | 140052 |
| Tested Applications | FIA |
| Alternative Names | Ethyl glucuronide, Ethanol |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid, may turn slightly opaque during storage |
| Concentration | 1.0 mg/ml (+/-10 %) (A280 nm, 1 mg/ml, 1 cm = 1.6) |
| Storage | 2-8 °C |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥ 95 % |
| Buffer | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, with 0.05 % sulfobetaine and 0.095 % NaN ₃ as a preservative |
| IEF Profile | 7.0-7.8 |
| Cross Reactivity | Not determined (N/D) |
| Specificity | Antibody recognizes ethyl glucuronide (EtG) |
| Matched Pairs | EtG-BSA antigen LA189 |
| Shelf Life | Unspecified |