

Product Datasheet

Anti-6-MAM R12703 SPTN-5 140041

Product Name	Anti-6-MAM R12703 SPTN-5
Catalog Number	140041
Tested Applications	FIA, Lateral Flow
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/-10 %)
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Not Determined
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.4-6.6
Cross Reactivity	Morphine ≤ 0.01 %
Specificity	Antibody recognizes 6-monoacetylmorphine (6-MAM)
Matched Pairs	6-MAM-BSA Antigen 170031
Shelf Life	Unspecified