

## Product Datasheet

### MDMA-BSA Antigen LA794

<b>Product Name</b>	MDMA-BSA Antigen
<b>Catalog Number</b>	LA794
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Short Term: ( ≤ 2 weeks) 2-8°C. Long Term -20 °C
<b>Source</b>	Synthetic
<b>Buffer</b>	20 mM Phosphate Buffered Saline, pH 7.2
<b>Matched Pairs</b>	HM869
<b>Shelf Life</b>	5 years