

## Product Datasheet

### Mycobacterium tuberculosis MPT64 (TB-MPT64) Antigen, Recombinant $\geq 95\%$ LA811

<b>Product Name</b>	Mycobacterium tuberculosis MPT64 (TB-MPT64) Antigen, Recombinant $\geq 95\%$
<b>Catalog Number</b>	LA811
<b>Tested Applications</b>	LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	$\leq -20^{\circ}\text{C}$
<b>Purity</b>	$\geq 95\%$ by SDS-PAGE
<b>Source</b>	E. coli
<b>Buffer</b>	Phosphate Buffered Saline, pH $7.0 \pm 0.2$
<b>Matched Pairs</b>	HM1423, HM1424
<b>Expression System</b>	E. coli