

## Product Datasheet

### **Mycobacterium tuberculosis MPT64 (TB-MPT64) Antibody HM1423**

<b>Product Name</b>	Mycobacterium tuberculosis MPT64 (TB-MPT64) Antibody
<b>Catalog Number</b>	HM1423
<b>Tested Applications</b>	LF
<b>Immunogen</b>	M. tuberculosis MPT64 (recombinant)
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥95% by SDS-PAGE
<b>Buffer</b>	Phosphate Buffered Saline, pH 7.0
<b>Matched Pairs</b>	HM1424 (detection)
<b>Preservatives</b>	0.02% Sodium Azide