

Product Datasheet

Brain Natriuretic Peptide (BNP) Antigen $\geq 90\%$ LA788

Product Name	Brain Natriuretic Peptide (BNP) Antigen $\geq 90\%$
Catalog Number	LA788
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Lyophilized Antigen
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -70°C. Avoid repeated freezing and thawing.
Purity	$\geq 90\%$
Source	Synthetic
Matched Pairs	HM818 and HM819