

Product Datasheet

Anti-Human CD62E Monoclonal Antibody, Azide Free Clone B-S3 853.810

Product Name	Anti-Human CD62E Monoclonal Antibody, Azide Free Clone B-S3
Catalog Number	853.810
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	Activated human endothelial cells
Alternative Names	ELAM-1 E-Selectin
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises the Endothelial Leucocyte Adhesion Molecule-1 (ELAM-1), E-selectin, a 97 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells