

# Product Datasheet

## Methamphetamine Antibody HM1268

Product Name	Methamphetamine Antibody
Catalog Number	HM1268
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	DEAE Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Amphetamine 0.8%, MDMA 0.3 %
Specificity	Methamphetamine
Matched Pairs	LA469 (capture)
Preservatives	0.05% Sodium Azide