

Product Datasheet

Clozapine-BSA Antigen LA800

| | |
|----------------------------|---|
| Product Name | Clozapine-BSA Antigen |
| Catalog Number | LA800 |
| Tested Applications | LF, ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | Liquid; Antigen-BSA conjugate |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.2 |
| Matched Pairs | HM1272 or HM1273 (detection) |
| Preservatives | 0.01% Sodium Azide |