

## Product Datasheet

### Clozapine-BSA Antigen LA801

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Clozapine-BSA Antigen   |
| <b>Catalog Number</b>      | LA801   |
| <b>Tested Applications</b> | LF, ELISA   |
| <b>Brand</b>               | EastCoast Bio   |
| <b>Form/Appearance</b>     | Liquid; Antigen-BSA conjugate   |
| <b>Storage</b>             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline, pH 7.2                                   |
| <b>Matched Pairs</b>       | HM1272 or HM1273 (detection)  |
| <b>Preservatives</b>       | 0.01% Sodium Azide  |