

Product Datasheet

Olanzapine Antibody HM1275

Product Name	Olanzapine Antibody
Catalog Number	HM1275
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Olanzapine-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90 %
Buffer	0.015M Potassium Phosphate Buffer with 0.85% NaCl, pH 7.2
Cross Reactivity	Olanzapine 100%; Aripiprazole, Clozapine, Olanzapine, Quetiapine, Risperidone, Ziprasidone, Amitriptyline, Clomipramine, Cyclobenzaprine, Desipramine, Doxepin, Imipramine, Maprotiline, Nordoxepin, Nortriptyline, Trimipramine <1%
Specificity	Olanzapine
Matched Pairs	LA802 (capture)
Preservatives	0.05% Sodium Azide