

# Product Datasheet

## Serum, Total Triiodothyronine (T3) Samples 991-24-PS-TT3

|                 |  |
|-----------------|--|
| Product Name    | Serum, Total Triiodothyronine (T3) Samples   |
| Catalog Number  | 991-24-PS-TT3  |
| Brand           | Lee Biosolutions   |
| Form/Appearance | Liquid   |
| Storage         | -20°C  |
| Source          | Single Human Donor   |
| Shelf Life      | Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time. |