

Product Datasheet

Serum, Triiodothyronine (T3) Uptake Samples 991-24-PS-UT3

| | |
|-----------------|--|
| Product Name | Serum, Triiodothyronine (T3) Uptake Samples |
| Catalog Number | 991-24-PS-UT3 |
| Brand | Lee Biosolutions |
| Form/Appearance | Liquid |
| Storage | -20°C |
| Source | Single Human Donor |
| Shelf Life | Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time. |