

Product Datasheet

Serum, Single Donor, Total Thyroxine (T4) 991-24-PS-TT4

Product Name	Serum, Single Donor, Total Thyroxine (T4)
Catalog Number	991-24-PS-TT4
Description	Each sample is a minimum of 1mL per donor. Inquire for larger volumes per donor or custom criteria.
Tested Applications	QC Materials
Alternative Names	Serum, fraction, clotted, post-coagulation, coagulated, powder, sera, blood serum, defibrinated, T4, Total T4, thyroxine
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	-20°C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None added, inquire if needed