

Product Datasheet

Serum, Troponin I Samples 991-24-PS-TNI

Product Name	Serum, Troponin I Samples
Catalog Number	991-24-PS-TNI
Species/Host	Human
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	-20 °C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.