

Product Datasheet

Anti-Human CD141 Monoclonal Antibody, Azide Free Clone B-S24 853.850.005

Product Name	Anti-Human CD141 Monoclonal Antibody, Azide Free Clone B-S24
Catalog Number	853.850.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	Eahy 926 cell lines
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Recombinant antibody from the original clone
Specificity	Recognises thrombomodulin (TM), a 105 kDa protein