

Product Datasheet

Anti-h PTH 10702 SPTN-5 101002

Product Name	Anti-h PTH 10702 SPTN-5
Catalog Number	101002
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Amino acids 1-2 of human PTH
Purity	≥ 95 %
Associated Products	Recombinant Parathyroid Hormone (PTH) Antigen 439-66
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.7-7.3
Specificity	Antibody recognizes human parathyroid hormone (PTH)
Matched Pairs	Anti-h PTH 10704, 10705, 10750
Shelf Life	Unspecified